BOSCH

vivalytic

SARS-CoV-2
*Severe Acute Respiratory Syndrome
Coronavirus 2*

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Package Contents

15 Vivalytic SARS-CoV-2 test cartridges for the detection of SARS-CoV-2 specific E gene sequence.

Safety Information

These instructions for use contain test-specific information only. For additional warnings and instructions refer to the instructions for use provided with your Vivalytic analyser (chapter device safety information). Only use Vivalytic cartridges and accessories approved for the analyser. Make sure to insert samples properly. Remove spilled samples with a cloth soaked with the appropriate solution. You can use 70 % ethanol. If you need to disinfect the cartridge by spraying (for example after heavy spillage of sample) do not use the cartridge afterwards since it might be damaged by disinfectant.

WARNING

- Do not use a cartridge if the sealed pouch or the cartridge itself is visibly damaged.
- Do not touch or scratch the detection area of the cartridge.
- Do not reuse a cartridge.
- Do not use expired cartridges. The expiration date can be found on the packaging and the cartridge label.
- Do not wait longer than 15 minutes after opening the cartridge pouch to begin the test. This maintains hygiene and avoids performance loss due to humidity. Prolonged exposure to humidity has a negative impact on test performance.
- Do not use sample types, media and volumes that are not approved for the test.
- Do not shake a cartridge that contains a sample.
- Do not turn the cartridge upside down.
- Handle potentially infectious patient samples and cartridges according to national laboratory standards and dispose samples and cartridges according to regional and laboratory standards.
Make sure to wear appropriate personal protective equipment (PPE).
- Place the cartridge on a clean surface.
- Be compliant with the national safety regulations and practices.

Additional Equipment required (required but not provided)

Bosch Vivalytic one Analyser Pipette (1000 µl)

Intended Use

The Bosch Vivalytic SARS-CoV-2 test is a qualitative PCR-based assay for the detection of SARS-CoV-2 from human nasopharyngeal and oropharyngeal swabs to aid in the diagnosis of a SARS-CoV-2 infection. Positive and negative results must be clinically correlated with patient history and other diagnostic information is necessary to determine patient infection status. Intended for use with a Vivalytic analyser by healthcare professionals.

Test Principle

Vivalytic SARS-CoV-2 is a real-time based PCR test utilizing reverse transcription (RT).

Storage und Usage Conditions

Product is stable until the expiry date if stored at +15 °C to 25 °C. Storage and usage conditions can be taken either from the cartridge, pouch or box label.

Reagents

All reagents necessary for the sample processing are integrated into the cartridge. The processing includes cell lysis, nucleic acid extraction, reverse transcription, DNA amplification and detection.

Reagents are PCR bead, binding buffer, washing buffer and elution buffer.

Sample Type/Medium

The test is intended for use with nasopharyngeal and oropharyngeal swab samples in eNAT™ medium (COPAN Italia S.P.A.). Collect and store samples as indicated in the manufacturer's data sheet.

Vivalytic SARS-CoV-2 – Instructions for use

Sample Preparation

Shake the sample tube containing the swab sample and the liquid medium for homogenization and fill 300 µl of homogenized patient sample in the sample input of the cartridge. Do not use viscous samples that are difficult to pipette.

Test Result

The test result (list of positively/negatively detected pathogens) is shown on the screen. For further information, you are guided through different levels of detail by the software. The PCR curve data level shows software-modified PCR curves. The curves are classified as positive or negative by the software. In case of positive curves, the respective C_q value is displayed.

Detection of the human control shows a successful extraction procedure and excludes an inhibition of the PCR reaction.

A test is rated as invalid if both target RNA and human control are not detected. Possible reasons for an invalid run might be poor sample quality due to a partial or complete absence of human cellular material in the sample.

Results are displayed for an invalid test but are not allowed to be used for diagnostic interpretation.

In the printed test report, pathogen, results, control and information on user, patient and analyser are listed with a signature field. In case of positive results, the respective C_q value is displayed on the report.

Pay attention to use the correct sample type, sample collection and storage of the sample and cartridges prior to the test run. Repeat the analysis with a new aliquot of the same sample if required.

In case of a failed test, first check for correct operating conditions of the analyser (refer to analyser's instructions for use, chapters device safety information and technical data). Restart the analyser. If the problem persists, contact the customer service.

As soon as a valid, positive result is shown on the screen, the user has the option to finish the test. Pending results are not displayed!

Results/Control

SARS-CoV-2	Human control	Validity	Result
+	+/-	valid	Sample is considered SARS-CoV-2 positive.
-	+	valid	Sample is considered negative.
-	-	invalid	Not evaluable. ¹

¹ Retesting is recommended.

Quality Control

If required by your local or laboratory standards, quality control testing has to be performed. You can either use pre-characterized patient samples that were investigated by a reference testing method or purchase quality control materials (e. g. AMPLIRUN® SARS-CoV-2 RNA CONTROL, Ref. MBC137-R, Vircell S. L. Spain). In case of unexpected results, repeat the analysis with another sample. If the result of a quality negative control sample e. g. pure eNAT™ medium remains positive, the analyser or its environment might be contaminated. Stop using the analyser and call the customer service. In case of repeated negative results for positive quality control samples, also call the customer service.

Limitations

The results of the Vivalytic SARS-CoV-2 test should be interpreted by a trained healthcare professional. The results of the Vivalytic SARS-CoV-2 test should not be used as the sole parameter for diagnosis.

- A negative result does not exclude pathogens being present in the sample at a level below assay sensitivity or a pathogen that is not covered by this assay.

- There is a risk of false negative values resulting from improperly collected, transported or handled samples.

Analytical Sensitivity (Limit of Detection, 95 % Detection rate)

The concentration at a detection rate of 95 % (LoD) was determined ([table 1](#)).

Inclusivity and Exclusivity

Specificity was ensured by the selection of primers and probes and their in silico analysis for possible cross-reactions based on publicly available nucleic acid sequences derived from the NCBI database.

To evaluate inclusivity, positive SARS-CoV-2 material ([table 2](#)) were purified and spiked into a PCR reaction and processed by a workflow using liquid components.

To exclude cross-reactivity (exclusivity), various strains of microorganisms representing common respiratory pathogens or closely related species were tested ([table 3](#)) by a workflow using liquid components. There was no evidence of microbial interference.

Interferences

Interferences were evaluated for endogenous and exogenous substances ([table 4](#)) that are potentially present in the patient sample. No interferences were detected.

Sensitivity and Specificity

The results derived from patient samples (positive and negative samples in Copan eNAT™ medium) collected in a clinical setting were compared with those of a reference method ([table 5](#)). In total, 69 samples were analysed.

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Table 4 – Tested Substances for Interference

Nasivin® spray (0.5 mg/ml oxymetazoline-hydrochloride); 5 % (v/v)

Rhinomer® (9 g/l NaCl); 10 % (v/v)

Weleda Heuschnupfenspray; 10 % (v/v)

Turixin® 2 % Mupirocin; 2.5 mg/ml

TobraZid® (80 mg Tobramycin / 2 ml ampoule); 4 µg/mL

Simulated nasal matrix [(2.5 % (w/v) porcine mucin, 1 % (v/v) human whole blood in 0.85 % sodium chloride (NaCl), in 1x PBS with 15 % glycerol]; 2.0 % (v/v)

Salbutamol-ratiopharm®; 5 % (v/v)

Neosynephrin-POS® 10 % eye drops; 5 % (v/v)

Otri-Allergie Spray Fluticason; 5 % (v/v)

Table 5 – Clinical Sensitivity [1] and Specificity [2] (95 % Confidence Interval)

SARS-CoV-2
[1] 97.6 % (87.4 – 99.9 %)
[2] 100 % (87.2 – 100.0 %)*

*(1) specimen was tested as “false positive” and was re-tested. Repetition of this specimen with a third reference assay yielded a positive result and the specimen was considered as “true positive” in this table. Without retesting specificity was 96.4 % (81.7 – 99.9 %).

Table 6 – Document History

Revision 01	Initial document
Revision 02	<u>Chapter - Test Result</u> <ol style="list-style-type: none">1. Introduction of PCR curve display description.2. Introduction of the possibility to finish a test earlier in case of positive results.3. Introduction of a tabular overview of results. <u>Annex – Performance Data</u> <ol style="list-style-type: none">1. Adaptation of Limit of Detection.2. Adaptation of concentration data in table 4 (simulated nasal matrix).

For more information see
www.bosch-vivalytic.com

Bosch Healthcare Solutions GmbH
Stuttgarter Straße 130
71332 Waiblingen, Germany

F09G 300 411

F09G300621_02